
## Top Glove, Top Quality, Top Efficiency Good Health, Safety First & Be Honest

## QA Department

Dumbbell Die Specification and Aging Method Guideline (01-00)

## a. Dumbbell Specification



Standard Reference	ASTM D3578	ASTM D6319	ASTM D3577	ASTM D4679	13392		ISO 11193/ KS ISO 11193 ISO 10282		MCAZ Zimbabwe		AS/NZS 4011		AS/NZS 4179		EN 455			
Glove Type	NR	NBR	NR SG	NR/NBR	NR	NBR	NR	NBR	NR SG	NR	NBR	NR	NBR	NR SG	NR	NBR	NR SG	
Dumb bell reference	ASTM D412						ISO 37							ASTM D412				
Dumb ben feference	Type C							Type II								Type D		
Overall length(mm) / a	min 115									min 100								
Overall width (mm) / b	25 <u>+</u> 1							12.5 <u>+</u> 1								16 <u>+</u> 1		
Narrow length (mm) / c	33 <u>+</u> 2									33 <u>+</u> 2								
Narrow width (mm) / d	6 <u>±</u> 0.05							4 <u>+</u> 0.1								3 <u>+</u> 0.05		
Test length (mm) / e	25 <u>+</u> 0.25							20 <u>+</u> 0.5								25 <u>+</u> 0.25		
Thickness (mm)	median average						median								average			

b. Aging Method

Standard Reference	ASTM D3578	ASTM D6319	ASTM D3577	ASTM D4679	NBR 13392 (ABNT)		ISO 11193		ISO 10282	MCAZ Zimbabwe			S/NZS 4011 AS/NZS 4179			EN 455:2	
Glove Type	NR	NBR	NR SG	NR/NBR	NR	NBR	NR	NBR	NR SG	NR	NBR	NR	NBR	NR SG	NR	NBR	NR SG
Aging Reference	ASTM D573									EN 455:2							
Aging Temperature / Period			0 <u>±</u> 2 °C /16 0 <u>±</u> 2 °C /22	66 <u>+</u> 2 Hrs 2 <u>+</u> 0.3 Hrs		70±2 °C / 168±2 hrs									70±2 °C / 168±2 hrs		
	1-Cut Sample									1-pack in inner							
Aging Method	2-Hanging Test Piece without touching each other							2-Ha	nging Tes	2-aging in oven							
			3-Conduc	ct Test	3-Conduct Test									3- conduct test			

c. Requirements

Standard Reference		ASTM D3578	ASTM D6319	ASTM D3577	ASTM D4679	NBR 13392 (ABNT)		ISO 11193/ KS ISO 11193		ISO 10282	MCAZ Zimbabwe		AS/NZS 4011		AS/NZS 4179			EN 155:2	
Glove Type		NR	NBR	NR SG	NR/NBR	NR	NBR	NR	NBR	NR SG	NR	NBR	NR	NBR	NR SG	NR	NBR	NR SG	
Tensile	Before Aging	18	14	24	10	21	12						21	12	23				
Strength (Mpa)	After Aging	14	14	18	7.5	16	9						16	-	17				
Force at	Before Aging							7	7	12.5	7	7				6	6	12	
Break (N)	After Aging							6	7	9.5	6	7				6	6	9	
Elongation (%)	Before Aging	650	500	750	500	700	500	650	500	700	650	500	700	300	700				
Eloligation (%)	After Aging	500	400	560	375	500	350	500	400	550	500	400	500	-	560				
Disposition		Acc: 1 Rej: 2												median					

Prepared by, Noor Azura Binti Azman/QA/MQC Asst Sup **Verified by,** Al-Fadilah Mohamad/QA Executive II Approved by, Noor Akilah Saidin/QA Deputy General Manager